


Anti-Kntc1 Rabbit Polyclonal Antibody For Wb, Elisa - P50748

Item Number: CM0014916



Introduction

Rabbit polyclonal antibody against human KNTC1, validated for WB and ELISA. Reacts with human and rat. Targets kinetochore-associated protein 1 essential for mitotic checkpoint. Suitable for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KNTC1 Rabbit pAb
Remarks/Alias	ROD; KNTC1
Species	Human
Gene ID	9735
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1100 of human KNTC1 (NP_055523.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P50748
Protein Weight	251kDa
Shipping	Ice bag