

Anti-[KO Validated] Calpain 2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17655

Item Number: CM0024023



Introduction

Rabbit polyclonal antibody targeting Calpain-2 catalytic subunit for WB, IHC-P, ELISA. Validated in human, mouse, rat. Suitable for cytoskeletal remodeling and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Calpain 2 Rabbit pAb
Remarks/Alias	CANP2, mCANP, CANPL2, CANPml, 2
Species	Human
Gene ID (Human)	824
Gene ID	824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-680 of human Calpain 2 (NP_001739.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17655
Protein Weight	80kDa
Shipping	Ice bag