

Anti-[Ko Validated] Edf1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Jmg1

Item Number: CM0020338

Introduction

Rabbit polyclonal antibody against mouse Edf1, knockout validated, suitable for Western blot and ELISA. Detects human and mouse Edf1, a transcriptional coactivator involved in endothelial differentiation and metabolism.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Edf1 Rabbit pAb
Remarks/Alias	MBF1; 0610008L11Rik; f1
Species	Mouse
GenelD (Human)	8721
GenelD	59022
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-148 of mouse Edf1 (NP_067494.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9JMG1
Protein Weight	16kDa
Shipping	Ice bag