

Anti-[Ko Validated] Pd-1/Cd279 Rabbit Polyclonal Antibody For Wb, Elisa - Q15116

Item Number: CM0029070



Introduction

Rabbit polyclonal antibody against human PD-1/CD279, knockout validated. Suitable for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying immune checkpoint pathways and cancer immunotherapy.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PD-1/CD279 Rabbit pAb
Remarks/Alias	PD1; PD-1; CD279; SLEB2; hPD-1; hPD-I; hSLE1; 79
Species	Human
GeneID (human)	5133
GeneID	5133
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15116
Protein Weight	32kDa
Shipping	Ice bag