

Anti-Kpna1 Rabbit Polyclonal Antibody For Wb And Elisa - P52294

Item Number: CM0014778



Introduction

Rabbit anti-KPNA1 polyclonal antibody for WB and ELISA; cross-reacts with human and mouse. Validated with recombinant immunogen; ideal for nuclear import research involving importin alpha-5.

[Learn More](#)

Parameter	Value
Product Name	KPNA1 Rabbit pAb
Remarks/Alias	RCH2; SRP1; IPOA5; NPI-1; KPNA1
Species	Human
GeneID (human)	3836
GeneID	3836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human KPNA1 (NP_002255.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52294
Protein Weight	60kDa
Shipping	Ice bag