

Anti-Kpna1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P52294

Item Number: CM0023511



Introduction

Rabbit polyclonal antibody against KPNA1 (importin subunit alpha-5). Validated for Western blot, IHC-P, and ELISA; reacts with human, mouse, and rat. Ideal for nuclear transport studies.

[Learn More](#)

Parameter	Value
Product Name	KPNA1 Rabbit pAb
Remarks/Alias	RCH2; SRP1; IPOA5; NPI-1; KPNA1
Species	Human
Gene ID (Human)	3836
Gene ID	3836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human KPNA1 (NP_002255.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52294
Protein Weight	60 kDa
Shipping	Ice bag