


Anti-Kpna2 Rabbit Monoclonal Antibody For Wb And Elisa - P52292

Item Number: CM0005113



Introduction

KPNA2 Rabbit Monoclonal Antibody validated for WB and ELISA with cross-reactivity to Human, Mouse, Rat. Recognizes Importin subunit alpha-1 (P52292), a nuclear transport adapter protein, ideal for studying nucleocytoplasmic trafficking.

[Learn More](#)

Parameter	Value
Product Name	KPNA2 Rabbit mAb
Remarks/Alias	QIP2; RCH1; IPOA1; PTAC58; SRP1alpha; SRP1-alpha; KPNA2
Species	Human
GeneID (Human)	3838
GeneID	3838
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KPNA2 (P52292).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52292
Protein Weight	58kDa
Shipping	Ice bag