

Anti-Kpna3 Rabbit Polyclonal Antibody For Wb And Elisa - 000505

Item Number: CM0019915

Introduction

Rabbit polyclonal antibody targeting human KPNA3 (importin alpha-4), suitable for Western blot and ELISA applications. Ideal for studying nuclear protein import, host-virus interactions, and hereditary spastic paraplegia.

[Learn More](#)

Parameter	Value
Product Name	KPNA3 Rabbit pAb
Remarks/Alias	SRP1; SRP4; IPOA4; SPG88; hSRP1; SRP1gamma; KPNA3
Species	Human
GeneID (Human)	3839
GeneID	3839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human KPNA3 (NP_002258.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00505
Protein Weight	58kDa
Shipping	Ice bag