

Anti-Kpna3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00505

Item Number: CM0030311

Introduction

High-quality KPNA3 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, and ELISA applications. Recognizes KPNA3 in human, mouse, rat. Targets importin subunit alpha-4, key in nuclear protein import and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	KPNA3 Rabbit pAb
Remarks/Alias	SRP1, SRP4, IPOA4, SPG88, hSRP1, SRP1gamma, KPNA3
Species	Human
Gene ID (Human)	3839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human KPNA3 (NP_002258.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00505
Protein Weight	58kDa
Shipping	Ice bag