

Anti-Kpna4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00629

Item Number: CM0024048



Introduction

Rabbit polyclonal antibody targeting human KPNA4 (Importin subunit alpha-3), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for nuclear transport and import assays.

[Learn More](#)

Parameter	Value
Product Name	KPNA4 Rabbit pAb
Remarks/Alias	QIP1; SRP3; IPOA3; KPNA4
Species	Human
Gene ID (Human)	3840
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human KPNA4 (NP_002259.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00629
Protein Weight	58kDa
Shipping	Ice bag