

Anti-Kpna5 Polyclonal Antibody For Wb, Elisa - 015131

Item Number: CM0029347



Introduction

KPNA5 Rabbit pAb is a polyclonal antibody raised in rabbit against human KPNA5 (importin subunit alpha-6). Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for nuclear transport studies.

[Learn More](#)

Parameter	Value
Product Name	KPNA5 Rabbit pAb
Remarks/Alias	SRP6; IPOA6; KPNA5
Species	Human
Gene ID	3841
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human KPNA5 (NP_002260.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15131
Protein Weight	61 kDa
Shipping	Ice bag