

Anti-Kpna6 Monoclonal Antibody For Wb, If/Icc, Elisa - O60684

Item Number: CM0006991



Introduction

Rabbit monoclonal antibody specific for human KPNA6 (importin alpha-7). Validated for WB, IF/ICC, and ELISA; cross-reacts with human and mouse. Detects a ~60 kDa band. Suitable for nuclear protein import and transport studies.

[Learn More](#)

Parameter	Value
Product Name	KPNA6 Rabbit mAb
Remarks/Alias	IPOA7; KPNA6
Species	Human
GeneID (Human)	23633
GeneID	23633
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60684
Protein Weight	60 kDa
Shipping	Ice bag