

Anti-Kpna6 Monoclonal Antibody For Wb - O60684

Item Number: CM0006846



Introduction

KPNA6 Rabbit PolymAb® is a monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. It reacts with human, mouse, rat KPNA6 (importin alpha-7), a nuclear import adaptor protein. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	KPNA6 Rabbit PolymAb®
Remarks/Alias	IPOA7; KPNA6
Species	Human
Gene ID (Human)	23633
Gene ID	23633
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 69-169 of human KPNA6 (NP_036448.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60684
Protein Weight	60 kDa
Shipping	Ice bag