

Anti-Kpna6 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O60684

Item Number: CM0006987



Introduction

KPNA6 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, and rat. Targets importin alpha-7, a key nuclear transport adapter, enabling reliable nucleocytoplasmic transport studies.

[Learn More](#)

Parameter	Value
Product Name	KPNA6 Rabbit mAb
Remarks/Alias	IPOA7; KPNA6
Species	Human
GeneID (human)	23633
GeneID	23633
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 69-169 of human KPNA6 (NP_036448.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60684
Protein Weight	60 kDa
Shipping	Ice bag