

Anti-Kpna6 Rabbit Polyclonal Antibody For Wb, Elisa - O60684

Item Number: CM0025467



Introduction

High-quality rabbit polyclonal antibody targeting human importin subunit alpha-7 (KPNA6). Validated for Western blot and ELISA, cross-reacts with human and mouse, ideal for nuclear import studies.

[Learn More](#)

Parameter	Value
Product Name	KPNA6 Rabbit pAb
Remarks/Alias	IPOA7; KPNA6
Species	Human
GeneID (Human)	23633
GeneID	23633
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 277-536 of human KPNA6 (NP_036448.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60684
Protein Weight	60kDa
Shipping	Ice bag