

## Anti-Kpnb1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14974

Item Number: CM0007710



### Introduction

This KPNB1 rabbit monoclonal antibody (Q14974) is validated for Western blot, immunohistochemistry, and ELISA, and cross-reacts with human, mouse, and rat. Ideal for studying importin beta-1 mediated nuclear import.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | KPNB1 Rabbit mAb                                                                                      |
| Remarks/Alias    | IMB1; IPO1; IPOB; Impnb; NTF97; KPNB1                                                                 |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 3837                                                                                                  |
| GeneID           | 3837                                                                                                  |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 300-400 of human KPNB1 (NP_002256.2) |
| Source           | Rabbit                                                                                                |
| Category         | Monoclonal Antibodies                                                                                 |
| Application      | WB, IHC-P, ELISA                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                     |
| SWISS            | Q14974                                                                                                |
| Protein Weight   | 97kDa                                                                                                 |
| Shipping         | Ice bag                                                                                               |