

Anti-Kras+Hras+Nras Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P01111/P01112/P01116

Item Number: CM0007717



Introduction

Recombinant rabbit monoclonal antibody targeting KRAS, HRAS, and NRAS. Suitable for WB, IF/ICC, IP, ELISA in human, mouse, rat samples. Ideal for RAS-MAPK signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KRAS+HRAS+NRAS Rabbit mAb
Remarks/Alias	C-K-RAS; CFC2; K-RAS2A; K-RAS2B; K-RAS4A; K-RAS4B; KI-RAS; KRAS1; KRAS2; NS; NS3; RALD; RASK2; c-Ki-ras2; KRAS+HRAS+NRAS
Species	Human
GeneID (Human)	4893/3845/3845
GeneID	4893/3265/3845
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Ras (P01111).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01111/P01112/P01116
Protein Weight	21kDa
Shipping	Ice bag