

Anti-Krbox4 Polyclonal Antibody For Wb, Elisa - Q5Juw0

Item Number: CM0017238

Introduction

KRBOX4 Rabbit pAb detects human, mouse, rat KRAB domain-containing protein 4 (ZNF673, KRABD4) by WB and ELISA. Validated with peptide immunogen, predicted molecular weight 20kDa, ideal for chromatin and transcriptional research.

[Learn More](#)

Parameter	Value
Product Name	KRBOX4 Rabbit pAb
Remarks/Alias	ZNF673; KRBOX4
Species	Human
GeneID (Human)	55634
GeneID	55634
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human KRBOX4 (NP_060246.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5JUW0
Protein Weight	20kDa
Shipping	Ice bag