


Anti-Kremen1 Polyclonal Antibody For Wb, Elisa - Q96Mu8

Item Number: CM0015123



Introduction

Rabbit polyclonal antibody targeting human KREMEN1, validated for WB and ELISA. Detects KREMEN1 in human, mouse, and rat samples. Suitable for Wnt/ β -catenin signaling, development, and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	KREMEN1 Rabbit pAb
Remarks/Alias	KRM1; ECTD13; KREMEN; KREMEN1
Species	Human
Gene ID (Human)	83999
Gene ID	83999
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 180-380 of human KREMEN1 (NP_001034659.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96MU8
Protein Weight	52kDa
Shipping	Ice bag