


Anti-Krr1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13601

Item Number: CM0026006



Introduction

KRR1 Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA. Targets the small subunit processome component involved in ribosome biogenesis. Validated in human, mouse, rat. Ideal for nucleolar and rRNA processing research.

[Learn More](#)

Parameter	Value
Product Name	KRR1 Rabbit pAb
Remarks/Alias	HRB2; RIP-1; KRR1
Species	Human
Gene ID (Human)	11103
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-381 of human KRR1 (NP_008974.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13601
Protein Weight	44kDa
Shipping	Ice bag