


Anti-Krt10 Rabbit Polyclonal Antibody For Wb And Elisa - P13645

Item Number: CM0024753



Introduction

Rabbit polyclonal antibody targeting human keratin 10 (KRT10/CK10). Validated for Western blot and ELISA. For studies of intermediate filaments, epidermal differentiation, and keratinocyte biology.

[Learn More](#)

Parameter	Value
Product Name	KRT10 Rabbit pAb
Remarks/Alias	BIE; EHK; IHL; K10; KPP; BCIE; CK10; KRT10
Species	Human
Gene ID (Human)	3858
Gene ID	3858
Immunogen	Synthetic peptide of human KRT10 (aa 350-460, NP_000412.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13645
Protein Weight	59kDa
Shipping	Ice bag