

Anti-Krt3 Polyclonal Antibody For Wb, If/Icc, Elisa - P12035

Item Number: CM0013853

Introduction

Rabbit polyclonal antibody targeting human KRT3, a corneal epithelium-specific type II keratin that pairs with KRT12. Validated for WB, IF/ICC, and ELISA. Ideal for studying intermediate filament networks and corneal dystrophies.

[Learn More](#)

Parameter	Value
Product Name	Anti-KRT3 Polyclonal Antibody
Remarks/Alias	K3; CK3; MECD2; KRT3
Species	Human
GeneID (Human)	3850
GeneID	3850
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-390 of human KRT3 (NP_476429.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P12035
Protein Weight	64kDa
Shipping	Ice bag