

Anti-Krt80 Rabbit Polyclonal Antibody For If/Icc, If-P, Elisa - Q6Kb66

Item Number: CM0027833

Introduction

High-quality rabbit polyclonal antibody targeting human keratin 80 (KRT80). Validated for IF/ICC, IF-P, and ELISA with cross-reactivity to mouse and rat. Explore keratin intermediate filament dynamics and tissue-specific expression.

[Learn More](#)

Parameter	Value
Product Name	KRT80 Rabbit pAb
Remarks/Alias	KB20; KRT80
Species	Human
GeneID (Human)	144501
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 242-452 of human KRT80 (NP_872313.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6KB66
Protein Weight	51kDa
Shipping	Ice bag