

Anti-Krt81 Rabbit Polyclonal Antibody For Wb, Elisa - Q14533

Item Number: CM0030747

Introduction

High-quality rabbit polyclonal antibody targeting human KRT81 (Keratin-81) for Western blot and ELISA applications. Ideal for hair keratin intermediate filament studies; cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	KRT81 Rabbit pAb
Remarks/Alias	HB1; K81; Hb-1; KRTHB1; MLN137; ghHkb1; hHAKB2-1; KRT81
Species	Human
Gene ID	3887
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human KRT81 (NP_002272.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14533
Protein Weight	55kDa
Shipping	Ice bag