

## Anti-Ksr1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ivt5

Item Number: CM0015755



### Introduction

High-quality rabbit polyclonal antibody against human KSR1, validated for Western blot, immunofluorescence, and ELISA. Recognizes human, mouse, and rat KSR1. Ideal for MAPK signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | KSR1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | KSR; RSU2; KSR1                                                                                                   |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 8844                                                                                                              |
| GeneID           | 8844                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 320-420 of human KSR1 (NP_055053.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q8IVT5                                                                                                            |
| Protein Weight   | 102kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |