

Anti-Ktn1 Polyclonal Antibody For Wb And Elisa - Q86Up2

Item Number: CM0013209

Introduction

Rabbit polyclonal antibody targeting human KTN1 (kinectin), validated for WB and ELISA. Cross-reacts with mouse and rat. Useful for investigating kinesin-driven vesicle motility, integrin adhesion, and ER membrane biology.

[Learn More](#)

Parameter	Value
Product Name	KTN1 Rabbit pAb
Remarks/Alias	CG1; KNT; MU-RMS-40.19; KTN1
Species	Human
GeneID (human)	3895
GeneID	3895
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KTN1 (NP_001072989.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UP2
Protein Weight	156kDa
Shipping	Ice bag