

Anti-Ktn1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86Up2

Item Number: CM0027026

Introduction

Rabbit polyclonal antibody targeting human KTN1 (kinectin, kinesin receptor), validated for WB, IHC-P, and ELISA. Detects the ~156 kDa ER membrane protein involved in vesicle motility and cell adhesion. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	KTN1 Rabbit pAb
Remarks/Alias	CG1; KNT; MU-RMS-40.19; KTN1
Species	Human
GeneID (Human)	3895
GeneID	3895
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human KTN1 (NP_001072989.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q86UP2
Protein Weight	156kDa
Shipping	Ice bag