

Anti-Ku70 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P12956

Item Number: CM0002079



Introduction

Rabbit monoclonal antibody against human Ku70 (XRCC6) for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Key tool for DNA double-strand break repair and NHEJ pathway research.

[Learn More](#)

Parameter	Value
Product Name	Ku70 Rabbit mAb
Remarks/Alias	ML8, KU70, TLAA, CTC75, CTCBF, G22P1, Ku70
Species	Human
Gene ID (Human)	2547
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 500-609 of human Ku70 (P12956)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12956
Protein Weight	70kDa
Shipping	Ice bag