

Anti-Ku70 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P12956

Item Number: CM0015544



Introduction

High-quality rabbit polyclonal antibody against human Ku70 (XRCC6). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for DNA repair research and assay development.

[Learn More](#)

Parameter	Details
Product Name	Ku70 Rabbit pAb
Remarks/Alias	ML8; KU70; TLAA; CTC75; CTCBF; G22P1; Ku70
Species	Human
GeneID (Human)	2547
GeneID	2547
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 408-609 of human Ku70 (NP_001460.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12956
Protein Weight	70kDa
Shipping	Ice bag