

Anti-Ku70 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P12956

Item Number: CM0026912



Introduction

High-quality rabbit polyclonal antibody targeting human Ku70 (XRCC6), validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for DNA damage repair and NHEJ pathway research.

[Learn More](#)

Parameter	Value
Product Name	Ku70 Rabbit pAb
Remarks/Alias	ML8; KU70; TLAA; CTC75; CTCBF; G22P1; Ku70
Species	Human
Gene ID	2547 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 501-600 of human Ku70 (NP_001460.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12956
Protein Weight	70kDa
Shipping	Ice bag