

Anti-Ku80 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P27641

Item Number: CM0022098

Introduction

Rabbit polyclonal antibody targeting mouse Ku80, validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human and mouse. Useful for DNA repair and telomere research.

[Learn More](#)

Parameter	Value
Product Name	Ku80 Rabbit pAb
Remarks/Alias	Ku80; Ku86; CTC85; CTCBF
Species	Mouse
GeneID (Human)	7520
GeneID	22596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 543-732 of mouse Ku80 (NP_033559.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P27641
Protein Weight	83kDa
Shipping	Ice bag