

Anti-Ku80 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P13010

Item Number: CM0028570



Introduction

Rabbit polyclonal antibody targeting human Ku80 (XRCC5). Validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Detects an 83kDa band. Suitable for DNA repair, NHEJ, and telomere research.

[Learn More](#)

Parameter	Value
Product Name	Ku80 Rabbit pAb
Remarks/Alias	KU80; KUB2; Ku86; NFIV; KARP1; KARP-1; Ku80
Species	Human
Gene ID (Human)	7520
Gene ID	7520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 543-732 of human Ku80 (NP_066964.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13010
Protein Weight	83kDa
Shipping	Ice bag