


Anti-Ku80 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P13010

Item Number: CM0002787



Introduction

Ku80 Rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, ELISA in human samples. Targets DNA repair protein Ku80 / XRCC5 for NHEJ and V(D)J recombination research. IVD raw material from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	Ku80 Rabbit mAb
Remarks/Alias	KU80; KUB2; Ku86; NFIV; KARP1; KARP-1; Ku80
Species	Human
GeneID	7520
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 633-732 of human Ku80 (P13010)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P13010
Protein Weight	83kDa
Shipping	Ice bag