

Anti-Kyat3 Rabbit Polyclonal Antibody (Kat3/Ccbl2) For Ihc-P, If/Icc, Elisa - Q6Yp21

Item Number: CM0028290

Introduction

KYAT3 rabbit polyclonal antibody (pAb) validated for IHC-P, IF/ICC, ELISA. Reacts with human and mouse. Targets kynurenine aminotransferase III, involved in tryptophan metabolism. Ideal for neuroscience and metabolic research.

[Learn More](#)

Parameter	Specification
Product Name	KYAT3 Rabbit pAb
Remarks/Alias	KAT3; CCBL2; KATIII
Species	Human
GeneID	56267
Immunogen	Synthetic peptide corresponding to amino acids 295-394 of human KYAT3 (NP_001008661.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6YP21
Protein Weight	51kDa
Shipping	Ice bag