

Anti-Kynu Monoclonal Antibody For Wb, Elisa - Q16719

Item Number: CM0004198



Introduction

This rabbit monoclonal antibody targets human kynureninase (KYNu), a key enzyme in tryptophan catabolism. It is suitable for Western blot and ELISA, with immunogen derived from recombinant protein. Cross-reactivity: human. Detects bands at 34 kDa and 52 kDa.

[Learn More](#)

Parameter	Value
Product Name	KYNu Rabbit mAb
Remarks/Alias	KYNu; KYNuU; kynureninase
Species	Human
GeneID (Human)	8942
GeneID	8942
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-465 of human KYNu (NP_003928.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16719
Protein Weight	34kDa/52kDa
Shipping	Ice bag