

Anti-Kynu Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16719

Item Number: CM0027914



Introduction

High-quality rabbit polyclonal antibody targeting human KYNU (kynureninase). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying tryptophan metabolism and neurobiology.

[Learn More](#)

Parameter	Value
Product Name	KYNU Rabbit pAb
Remarks/Alias	KYNUU; VCRL2; KYNU
Species	Human
GeneID (Human)	8942
GeneID	8942
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-307 of human KYNU (NP_001028170.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16719
Protein Weight	52kDa
Shipping	Ice bag