

Anti-L-Myc Rabbit Polyclonal Antibody For Wb, Elisa - P12524

Item Number: CM0020020

Introduction

Rabbit polyclonal antibody against human L-Myc (MYCL) validated for western blot and ELISA. Suitable for studying transcription factor function in nuclear processes and oncogenic signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	L-Myc Rabbit pAb
Remarks/Alias	LMYC; L-Myc; MYCL1; bHLHe38; MYCL
Species	Human
GeneID (human)	4610
GeneID	4610
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human MYCL (NP_001028253.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P12524
Protein Weight	40kDa
Shipping	Ice bag