


Anti-L-Plastin/Lcp1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P13796

Item Number: CM0011135



Introduction

Rabbit monoclonal antibody against human L-Plastin/LCP1 (Plastin-2), validated for WB, IHC-P, ELISA. Suitable for studies of T-cell activation, actin cytoskeleton, and immune cell function.

[Learn More](#)

Parameter	Value
Product Name	L-Plastin/LCP1 Rabbit mAb
Remarks/Alias	LPL; CP64; PLS2; LC64P; HEL-S-37; L-PLASTIN
Species	Human
Gene ID (Human)	3936
Gene ID	3936
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 528-627 of human L-Plastin/LCP1 (P13796).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P13796
Protein Weight	70kDa
Shipping	Ice bag