


Anti-L-Plastin/Lcp1 Polyclonal Antibody For Wb, Ip, Elisa - P13796

Item Number: CM0030189



Introduction

Rabbit polyclonal antibody against human L-Plastin/LCP1 (Plastin-2, LCP-1). Validated for WB, IP, ELISA. Cross-reacts with human and mouse. Ideal for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	L-Plastin/LCP1 Rabbit pAb
Remarks/Alias	LPL; CP64; PLS2; LC64P; HEL-S-37; L-PLASTIN; L-Plastin/LCP1
Species	Human
GeneID	3936
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human L-Plastin/LCP1 (NP_002289.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13796
Protein Weight	70kDa
Shipping	Ice bag