

Anti-L1R Rabbit Polyclonal Antibody For Wb, Elisa - Q8V4Z7

Item Number: CM0027628

Introduction

Rabbit polyclonal antibody against Monkeypox virus L1R protein (strain Zaire-96-I-16).
Validated for WB and ELISA, cross-reactive with the immunizing strain. Suitable for monkeypox research and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	L1R Rabbit pAb
Species	Monkeypox virus Zaire-96-I-16
GeneID	928938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of Monkeypox virus Zaire-96-I-16 L1R (NP_536512.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Monkeypox virus Zaire-96-I-16
SWISS	Q8V4Z7
Protein Weight	18kDa
Shipping	Ice bag