

Anti-L2Hgdh Polyclonal Antibody For Wb, Elisa - Q9H9P8

Item Number: CM0027231

Introduction

Rabbit polyclonal antibody targeting human L2HGDH (L-2-hydroxyglutarate dehydrogenase), validated for WB and ELISA. Reacts with human and mouse. Ideal for mitochondrial metabolism and disease research.

[Learn More](#)

Parameter	Value
Product Name	L2HGDH Rabbit pAb
Remarks/Alias	L2HGA; C14orf160; L2HGDH
Species	Human
Gene ID	79944
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human L2HGDH (NP_079160.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H9P8
Protein Weight	50kDa
Shipping	Ice bag