


Anti-L3MbtL1 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q9Y468

Item Number: CM0029729



Introduction

This anti-L3MBTL1 rabbit polyclonal antibody is validated for IF/ICC and ELISA, targeting human L3MBTL1, a Polycomb group chromatin regulator that represses transcription via histone methylation recognition.

[Learn More](#)

Parameter	Value
Product Name	L3MBTL1 Rabbit pAb
Remarks/Alias	L3MBTL; ZC2HC3; H-L(3)MBT; dj138B7.3; L3MBTL1
Species	Human
GeneID (Human)	26013
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 533-772 of human L3MBTL1 (NP_056293.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y468
Protein Weight	92kDa
Shipping	Ice bag