

Anti-L3MbtI2 Polyclonal Antibody For Wb, Elisa - Q969R5

Item Number: CM0015614



Introduction

L3MBTL2 rabbit polyclonal antibody validated for Western blot and ELISA. Targets human L3MBTL2, a Polycomb group chromatin regulator with a molecular weight of 79 kDa.

[Learn More](#)

Parameter	Value
Product Name	L3MBTL2 Rabbit pAb
Remarks/Alias	L3MBT; H-I(3)mbt-I; L3MBTL2
Species	Human
GeneID	83746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 486-705 of human L3MBTL2 (NP_113676.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q969R5
Protein Weight	79kDa
Shipping	Ice bag