


Anti-L3MbtI3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Jm7

Item Number: CM0030927



Introduction

Rabbit polyclonal antibody against L3MBTL3 (MBT-1), validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse L3MBTL3. Ideal for Notch signaling, chromatin regulation, and ubiquitin ligase studies.

[Learn More](#)

Parameter	Value
Product Name	L3MBTL3 Rabbit pAb
Remarks/Alias	MBT1; MBT-1; L3MBTL3
Species	Human
Gene ID (Human)	84456
Gene ID	84456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human L3MBTL3 (NP_115814.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96JM7
Protein Weight	88kDa
Shipping	Ice bag