


Anti-Lactb Rabbit Polyclonal Antibody For Wb, Elisa - P83111

Item Number: CM0023935



Introduction

Anti-LACTB Rabbit Polyclonal Antibody (P83111) validated for WB and ELISA with human, mouse, rat reactivity. Targets mitochondrial serine protease and tumor suppressor LACTB; ideal for cancer and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	LACTB Rabbit pAb
Remarks/Alias	LACTB; G24; MRPL56; lactamase beta
Species	Human
Gene ID (Human)	114294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 116-255 of human LACTB (NP_116246.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P83111
Protein Weight	41kDa/60kDa
Shipping	Ice bag