


Anti-Lactyl-Histone H4 (K5) Rabbit Polyclonal Antibody (P62805) For Wb, Db, Elisa

Item Number: CM0000074



Introduction

Rabbit polyclonal antibody targeting lactylated histone H4 at lysine 5 (K5-lactylation). Suitable for Western blot, dot blot, and ELISA applications in human, mouse, and rat samples with wide cross-reactivity predicted. UniProt P62805.

[Learn More](#)

Parameter	Value
Product Name	Lactic acid -Histone H4-K5 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Lactic acid -Histone H4-K5
Species	Human
GeneID (human)	8370
GeneID	8359
Immunogen	A synthetic lactylated peptide around K5 of human Histone H4 (NP_003539.1).
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag