

Anti-Ladinin-1 (Lad1) Rabbit Polyclonal Antibody For Wb And Elisa - O00515

Item Number: CM0021690

Introduction

Validated LAD1 rabbit polyclonal antibody for Western blotting and ELISA, with reactivity against human, mouse, and rat samples. Ideal for basement membrane and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	LAD1 Rabbit pAb
Remarks/Alias	LadA; LAD1
Species	Human
GeneID (human)	3898
GeneID	3898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-200 of human LAD1 (NP_005549.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00515
Protein Weight	57kDa
Shipping	Ice bag