


Anti-Lag3 Rabbit Polyclonal Antibody For Wb And Elisa - P18627

Item Number: CM0025217



Introduction

Rabbit polyclonal antibody targeting human LAG3 (CD223), validated for Western blot and ELISA. Cross-reacts with mouse and rat; suitable for immune checkpoint studies and cancer immunotherapy research.

[Learn More](#)

Parameter	Value
Product Name	LAG3 Rabbit pAb
Remarks/Alias	CD223; LAG3
Species	Human
GeneID (Human)	3902
GeneID	3902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-170 of human LAG3 (NP_002277.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18627
Protein Weight	57kDa
Shipping	Ice bag