

Anti-Lair1 Rabbit Polyclonal Antibody For Wb And Elisa - Q6Gtx8

Item Number: CM0017448



Introduction

Rabbit polyclonal antibody against human LAIR1 (CD305), validated for WB and ELISA. Targets an inhibitory immune receptor on NK, B, and T cells. Suitable for immunology and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	LAIR1 Rabbit pAb
Remarks/Alias	CD305; LAIR-1; LAIR1
Species	Human
Gene ID	3903
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-163 of human LAIR1 (NP_002278.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6GTX8
Protein Weight	31kDa
Shipping	Ice bag