


Anti-Lalba Polyclonal Antibody For Wb, Ihc-P, Elisa - P00709

Item Number: CM0027115



Introduction

Rabbit polyclonal antibody targeting human LALBA (alpha-lactalbumin). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Useful for lactose synthase studies and mammary gland biology.

[Learn More](#)

Parameter	Value
Product Name	LALBA Rabbit pAb
Remarks/Alias	LYZG; HAMLET; LALBA
Species	Human
GeneID (Human)	3906
GeneID	3906
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-142 of human LALBA (NP_002280.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P00709
Protein Weight	16kDa
Shipping	Ice bag